



Avalus™

PERIGON Pivotal Trial 5-year clinical update

The PERIGON Pivotal Trial was designed to evaluate the safety and effectiveness of the Model 400 (Avalus) bovine pericardial stented aortic bioprosthesis in a patient population undergoing surgical aortic valve replacement. Patients were from Canada (268), Europe (475) and the United States (375). Of the 1,288 patients enrolled, 1,118 were implanted with the Avalus valve. At 5 years, the Avalus valve has shown 0 SVD.

Key Takeaways

The PERIGON Pivotal Trial 5-year clinical update demonstrated safety and effectiveness of Avalus™ Bioprosthesis:


- EXCELLENT DURABILITY through 5 years: No SVD
- STABLE HEMODYNAMICS: Improvements in mean aortic gradients were stable over time in all valve sizes
- EXCELLENT FUNCTIONALITY: > 95% of subjects had less than mild central regurgitation at 5 years
- STABLE NYHA CLASS: > 95% of subjects reported NYHA class I or II at 5 years
- HIGH survival rate through 5 years

Baseline demographics

Patient characteristics	N = 1,118
Age, years	70.2 ± 9.0
Male sex	75.1%
BSA, m ²	2.0 ± 0.2
STS risk of mortality	2.0 ± 1.4%
NYHA class I/II	57.8%
NYHA class III/IV	42.2%
Atrial fibrillation	10.5%
Coronary artery disease	43.6%

Kaplan-Meier survival analysis

Freedom-from-event rate for all-cause mortality and structural valve deterioration




Procedural data

Procedural characteristics	N = 1,118
Primary indication	
Aortic stenosis	84.3%
Aortic regurgitation	5.7%
Mixed AS/AR	9.5%
Surgical approach	
Median sternotomy	79.6%
Less invasive approach	20.4%
Concomitant CABG	32.4%


Echocardiographic findings

All valve sizes†



†Data are not paired.


Valve size distribution



*Not commercially available in the U.S.

Echocardiographic findings

Mean aortic gradient, By valve size




STABLE LOW Mean Aortic Gradients in all valve sizes

N represents number of study patients implanted with mentioned valve size. Number of echocardiography images analyzed to date differ for each study timepoint.

Echocardiographic findings

Transvalvular regurgitation

All valve sizes†




> 95% of subjects had less than mild central regurgitation at 5 years

NYHA classification by visit

NYHA classification

All valve sizes†



> 95% of subjects reported NYHA class I or II at 5 years

Medtronic

Europe
Medtronic International Trading Sàrl.
Route du Molliau 31
Case postale
CH-1131 Tolochenaz
Tel: +41 (0)21 802 70 00
Fax: +41 (0)21 802 79 00

medtronic.eu

PERIGON Pivotal Trial data on file as of October 2021.

See the device manual for detailed information regarding the instructions for use, indications, contraindications, warnings, precautions, and potential adverse events. For further information, contact your local Medtronic representative and/or consult the Medtronic website at medtronic.eu.

For applicable products, consult instructions for use on manuals.medtronic.com. Manuals can be viewed using a current version of any major internet browser. For best results, use Adobe Acrobat® Reader with the browser.

UC202208173EE © Medtronic 2021. All rights reserved.