



# Endurant™ II/IIs

AAA Stent Graft System

## PRODUCT CODES

### ENDURANT II BIFURCATIONS

Product Code						
	Proximal Graft Diameter (mm)	Distal Graft Diameter (mm)	Distal Design	Total Covered Length (mm)	Delivery System	Catheter Outer Diameter (F)
ETBF	23	13	C	124	EE	18
ETBF	23	13	C	145	EE	18
ETBF	23	13	C	166	EE	18
ETBF	23	16	C	124	EE	18
ETBF	23	16	C	145	EE	18
ETBF	23	16	C	166	EE	18
ETBF	25	13	C	124	EE	18
ETBF	25	13	C	145	EE	18
ETBF	25	13	C	166	EE	18
ETBF	25	16	C	124	EE	18
ETBF	25	16	C	145	EE	18
ETBF	25	16	C	166	EE	18
ETBF	28	13	C	124	EE	18
ETBF	28	13	C	145	EE	18
ETBF	28	13	C	166	EE	18
ETBF	28	16	C	124	EE	18
ETBF	28	16	C	145	EE	18
ETBF	28	16	C	166	EE	18
ETBF	28	20	C	124	EE	18
ETBF	28	20	C	145	EE	18
ETBF	28	20	C	166	EE	18
ETBF	32	16	C	124	EE	20
ETBF	32	16	C	145	EE	20
ETBF	32	16	C	166	EE	20
ETBF	32	20	C	124	EE	20
ETBF	32	20	C	145	EE	20
ETBF	32	20	C	166	EE	20
ETBF	36	16	C	145	EE	20
ETBF	36	16	C	166	EE	20
ETBF	36	20	C	145	EE	20
ETBF	36	20	C	166	EE	20

### ENDURANT IIs BIFURCATIONS

Product Code						
	Proximal Graft Diameter (mm)	Distal Graft Diameter (mm)	Distal Design	Total Covered Length (mm)	Delivery System	Catheter Outer Diameter (F)
ESBF	23	14	C	103	EE	18
ESBF	25	14	C	103	EE	18
ESBF	28	14	C	103	EE	18
ESBF	32	14	C	103	EE	20
ESBF	36	14	C	103	EE	20

### LIMBS\*

Product Code								
	Proximal Graft Diameter (mm)	Distal Graft Diameter (mm)	Distal Design	Total Covered Length (mm)	Delivery System	Catheter Outer Diameter (F)	Total Contralateral Covered Length with Ell/EIIs Bifurcated <sup>†</sup>	Total Ipsilateral Covered Length with Ells Bifurcated <sup>†</sup>
ETLW	16	10	C	82	EE	14	136	155
ETLW	16	10	C	93	EE	14	147	166
ETLW	16	10	C	124	EE	14	178	177-197
ETLW	16	10	C	156	EE	16	210	209-229
ETLW	16	10	C	199	EE	16	253	252-272
ETLW	16	13	C	82	EE	14	136	155
ETLW	16	13	C	93	EE	14	147	166
ETLW	16	13	C	124	EE	14	178	177-197
ETLW	16	13	C	156	EE	16	210	209-229
ETLW	16	13	C	199	EE	16	253	252-272
ETLW	16	16	C	82	EE	14	136	135-155
ETLW	16	16	C	93	EE	14	147	146-166
ETLW	16	16	C	124	EE	14	178	177-197
ETLW	16	16	C	156	EE	16	210	209-229
ETLW	16	16	C	199	EE	16	253	252-272
ETLW	16	20	C	82	EE	16	136	155
ETLW	16	20	C	93	EE	16	147	166
ETLW	16	20	C	124	EE	16	178	177-197
ETLW	16	20	C	156	EE	16	210	209-229
ETLW	16	20	C	199	EE	16	253	252-272
ETLW	16	24	C	82	EE	16	136	155
ETLW	16	24	C	93	EE	16	147	166
ETLW	16	24	C	124	EE	16	178	177-197
ETLW	16	24	C	156	EE	16	210	209-229
ETLW	16	24	C	199	EE	16	253	252-272
ETLW	16	28	C	82	EE	16	136	155
ETLW	16	28	C	93	EE	16	147	166
ETLW	16	28	C	124	EE	16	178	177-197
ETLW	16	28	C	156	EE	16	210	209-229
ETLW	16	28	C	199	EE	16	253	252-272

\*The limb mates with the AUI stent graft on the ipsilateral side.

<sup>†</sup>These calculations assume the minimum 30 mm overlap between the bifurcated stent graft and the contralateral iliac limb per the Endurant II Stent Graft System Instructions For Use. When using the 124 mm length bifurcated stent graft, subtract 10 mm from Total Contralateral Covered Length with Bifurcated.

<sup>†</sup>The 3-5 stent overlap is available only with select limbs. Please refer to the Instructions For Use for more information.

### AORTIC EXTENSIONS

Product Code						
	Proximal Graft Diameter (mm)	Distal Graft Diameter (mm)	Distal Design	Total Covered Length (mm)	Delivery System	Catheter Outer Diameter (F)
ETCF	23	23	C	49	EE	18
ETCF	25	25	C	49	EE	18
ETCF	28	28	C	49	EE	18
ETCF	32	32	C	49	EE	20
ETCF	36	36	C	49	EE	20

### ABDOMINAL TUBES

Product Code						
	Proximal Graft Diameter (mm)	Distal Graft Diameter (mm)	Distal Design	Total Covered Length (mm)	Delivery System	Catheter Outer Diameter (F)
ETTF	23	23	C	70	EE	18
ETTF	25	25	C	70	EE	18
ETTF	28	28	C	70	EE	18
ETTF	32	32	C	70	EE	20
ETTF	36	36	C	70	EE	20

### AUI

Product Code						
	Proximal Graft Diameter (mm)	Distal Graft Diameter (mm)	Distal Design	Total Covered Length (mm)	Delivery System	Catheter Outer Diameter (F)
ETUF	23	14	C	102	EE	18
ETUF	25	14	C	102	EE	18
ETUF	28	14	C	102	EE	18
ETUF	32	14	C	102	EE	20
ETUF	36	14	C	102	EE	20

### ILIAC EXTENSIONS

Product Code						
	Proximal Graft Diameter (mm)	Distal Graft Diameter (mm)	Distal Design	Total Covered Length (mm)	Delivery System	Catheter Outer Diameter (F)
ETEW	10	10	C	82	EE	14
ETEW	13	13	C	82	EE	14
ETEW	20	20	C	82	EE	16
ETEW	24	24	C	82	EE	16
ETEW	28	28	C	82	EE	18

# Medtronic

**Invatec S.p.A.**

Via Martiri della Libertà 7 25030  
Roncadelle (BS)  
Italy  
Tel: +39.030.2589311

**Medtronic**

710 Medtronic Parkway NE  
Minneapolis, MN 55432  
USA  
Tel: +1.763.514.4000

**Medtronic International Trading  
Sàrl**

Route du Molliau 31  
CH-1131 Tolochenaz  
Switzerland  
Tel: +41.21.802.7000

**Medtronic of Canada Ltd.**

99 Hereford Street  
Brampton, Ontario L6Y 0R3  
Canada  
Tel: +1.905.460.3800

**Medtronic Latin America**

3750 NW 87 Avenue, Suite 700  
Miami, FL 33178  
USA  
Tel: +1.305.500.9328

**Medtronic International Ltd.**

49 Changi South  
Avenue 2  
Singapore 486056  
Tel: +65.6436.5000

**Medtronic Australasia Pty Ltd.**

97 Waterloo Road  
North Ryde, NSW 2113  
Australia  
Tel: +61.29857.9000

**Medtronic New Zealand**

Unit N16, Mezzanine Level 5  
Gloucester Park Road  
Onehunga  
Auckland  
New Zealand

**Medtronic Korea Co., Ltd.**

5F, Sajo Building  
1001 Daechi-dong, Kangnam-ku  
Seoul, 135-280  
Korea  
Tel: +82.2.3404.3600

**Medtronic META FZ-LLC**

Office Park, Block D, 2nd Floor  
PO Box 500638  
Dubai Internet City | Dubai  
United Arab Emirates  
Tel: +971.4.818.2666

[medtronic.com/aortic](http://medtronic.com/aortic)